

Suggested Work Flow & Reagents



1 SAMPLE COLLECTION



2 FIXATION

95% Ethanol Solution (catalog# 3340-G)



3 RINSE

Deionized water (catalog# 3315-G)



4 STAINING

Hematoxylin, Gill 1X (catalog# 7012-G)
Deionized water (catalog# 3315-G)
Bluing Reagent (catalog# 3356-G)
95% Ethanol Solution (catalog# 3340-G)
Pap Stain, OG-6 (catalog# 7042-G)
Pap Stain, EA-50 (catalog# 7028-G) or
Pap Stain, EA-36 (catalog# 7026-G)
100% Ethanol Solution (catalog# 3341-G)
Xylene (catalog# 3346-G)



5 MICROSCOPIC ANALYSIS

Acrylic Mounting Medium (Catalog# 3349-04,16)

Go to ASTRALDIAGNOSTICS.COM to see our full product listings.
Please contact us regarding special orders or to learn about our OEM capabilities.