



Our hemoglobin products serve as biologic quality controls for analyzer and assay development and laboratory-based QC controls (e.g. linearity controls). They have additionally been utilized as immunogens for antibody development.

### Hemoglobin Products

Cat #	Product	Amount
HbA0	Human Hemoglobin A0 (HbA0)	5.0 mg
HbA2	Human Hemoglobin A1a, A1b mixture	1.0 mg
HbA1	Human Hemoglobin A1c	1.0 mg
gly HbA0	Glycated Human Hemoglobin A0 (HbA0)	1.0 mg
de HbA0	Deglycated Human Hemoglobin A0 (HbA0)	5.0 mg

### Hemoglobin Assay

Cat #	Assay	Target/Species	Assay Principle	Specimen	Assay Range	# Plates
1044	Drabkin's Microplate Hemoglobin Assay	Hemoglobin (all species)	Chemical reaction, absorbance	Blood Samples, column fractions, etc.	0.5mg/ml - 20 mg/ml	1 x 96-well plate

Over 80 years of Experience Manufacturing Accurate, Repeatable, & Easy-to-use Diagnostics & Reagents

Let's Connect!



service@ethosbiosciences.com  
856.224.0900

