



Ethos Biosciences' Nephrology ELISAs measure proteins and biomarkers that are of significance for kidney disease and other metabolic disorders. They are suitable for use with a range of sample types, including urine, plasma, serum and cell culture supernatants. Kits feature a 96-well plate format, intended for up to 40 unique samples per plate (in addition to duplicates and controls).

The Creatinine Companion (Cat. 1012) is commonly used in conjunction with albumin ELISAs for additional intra- and inter-assay standardization.

Cat #	Assay	Target/Species	Assay Principle	Specimen	Assay Range	# Plates
1004	Albuwell Hu	Albumin [Human]	direct competitive	urine [Hu]	0.156 - 10µg/ml	2 x 96-well plates
1011	Albuwell M	Albumin [Mouse]	direct competitive	urine [Mouse]	0.156 - 10µg/ml	2 x 96-well plates
NR002	Nephurat II	Albumin [Rat]	direct competitive	urine [Rat]	2.0 - 100µg/ml	2 x 96-well plates
1012	Creatinine Companion	Creatinine [Human, Mouse, Rat]	Jaffe reaction	urine [all species]	10 - 100µg/ml	2 x 96-well plates
1035	Human Nephryn	Nephryn [Human]	indirect competitive	urine [Hu]	0.031 - 2.0µg/ml	1 x 96-well plate
1036	Rat Nephryn	Nephryn [Rat]	indirect competitive	urine [Rat]	0.156 - 10µg/ml	1 x 96-well plate
1037	Mouse Nephryn	Nephryn [Mouse]	indirect competitive	urine [Mouse]	0.031 - 2.0µg/ml	1 x 96-well plate
1041	Cystatin C, Human	Cystatin C [Human]	direct competitive	urine, serum, saliva, cell culture supernatant	7.8 - 500ng/ml	1 x 96-well plate
1042	Cystatin C, Rat	Cystatin C [Rat]	direct competitive	urine, serum, saliva, cell culture supernatant	7.8 - 500ng/ml	1 x 96-well plate
1043	Cystatin C, Mouse	Cystatin C [Mouse]	direct competitive	urine, serum, saliva, cell culture supernatant	7.8 - 500ng/ml	1 x 96-well plate
1038	Human NGAL	NGAL [Human]	direct [Ab] sandwich	urine, serum, cell culture [Hu]	0.2 - 12.5ng/ml	1 x 96-well plate
1039	Rat NGAL	NGAL [Rat]	direct [Ab] sandwich	urine, serum, cell culture [Rat]	0.4 - 25ng/ml	1 x 96-well plate
1040	Mouse NGAL	NGAL [Mouse]	indirect competitive	urine, serum, cell culture [Mouse]	3.12 - 200ng/ml	1 x 96-well plate

Over 80 years of Experience Manufacturing Accurate, Repeatable, & Easy-to-use Diagnostics & Reagents

Let's Connect!



service@ethosbiosciences.com  
856.224.0900

